KD-PENOFINE

Pen Needles

Pen-Nadeln • Pennkanyler • Pennåle • Pennekanyler • Insuliinikynäneulat • Pennaalden • Aiguilles pour stylos • Igły do pena • Jehly pro injekční pera • Ihly na perá • Igle za inzulinsku olovku • Ace stilou injector • Игли за писалки

Compatible pens


Kompatible Pens • Kompatibla pennor • Kompatible penne • Kompatible penner • Yhteensopivat kynät • Compatibele pennen • Stylos compatibles • Kompatybilne peny • Kompatibilní pera • Kompatibilné perá • Kompatibilne inzulinske olovke • Ace stilou injector compatibile • Съвместими писалки

BD Medical	Puregon Pen®
Berlin Chemie	BerliPen® areo 2 BerliPen® areo 3
Bioton	GensuPen Improve GensuPen2
Eli Lilly	Humalog® KwikPen™ HumaPen® HumaPen® Ergo HumaPen® Ergo II HumaPen® Luxura HumaPen® Luxura HD Saline KwikPen™ HumaPen® Savvio
Gan & Lee	GanleePen™
Novo Nordisk	NovoMix®30 FlexPen® NovoPen® 3 NovoPen® 4 NovoPen® 5 NovoPen® 6 NovoPen® Junior NovoPen Echo® NovoPen Echo® Plus Tresiba® FlexTouch™
Owen Mumford	AutoPen® AutoPen® 2 AutoPen® 24
Sanofi Aventis	ALLSTAR® ALLSTAR® Pro JuniorSTAR® Lantus® SoloStar® OptiPen® Pro 1 Sanofi Aventis Lantus® TactiPen®
Ypsomed	YpsoPen®

EN ISO 11608-2

Figures 1-9

Abbildung 1 bis 9 • Figur 1-9 • Figur 1-9 • Figur 1-9 • Kuvat 1-9 • Afbeelding 1-9 • Figures 1-9 • Rysunki 1-9 • Obrázky 1-9 • Obrázky 1 – 9 • Slike 1-9 • Figurile 1-9 • Фигури 1 – 9

en Pen Needles

Brief description:

KD-*PENOFINE* pen needles consist of the following components: big protective cap, sealing strip, small protective cap and the needle fixed in a plastic hub.

Materials used:

Stainless steel, polypropylene (PP), silicone oil (lubricant), UV glue (adhesive)

Intended purpose:

KD-*PENOFINE* pen needles are single use medical device intended to be connected to a compatible insulin pen for subcutaneous insulin administration within the frame of a diabetes mellitus treatment regimen. KD-*PENOFINE* pen needles are sterilized using ethylene oxide.

Indication:

Subcutaneous injections of insulin with a compatible insulin pen

Contraindications:

- Insulin/insulin derivatives injections during episodes of hypoglycemia
- Hypersensitivity to insulin/insulin derivatives or one of its excipients
- Inflammation, tumors, ulceration or other pathologic skin / subcutaneous tissue alterations or lesions such as bruising or indurations or scar tissue at the injection site

Residual risks, undesirable side effects:

- Injection-related side effects such as bleeding, bruising (hematoma) or skin reactions (e.g. redness)
- Lipohypertrophy
- Hypoglycemia
- Needle breaking or bending
- Pain
- Insulin backflow (leakage)

Intended user:

The use of this product is restricted to qualified healthcare professionals and trained laypersons such as patients or caregivers.

Intended patient population:

KD-*PENOFINE* pen needles are intended for patients who need insulin injection for blood glucose control without gender or age constraints. The different needle sizes available are to be selected by experienced healthcare professionals regarding suitability for the intended injection sites and the patient's body condition.

⚠ Warnings, precautions:

- Read instructions before use. KD-*PENOFINE* is to be used according to these instructions for use.
- **As a layperson, consult your doctor regarding the most suitable needle size, injection technique and injection sites to be used and always follow your doctor's instructions.**
- Check the compatibility between the pen needle and the insulin pen — a list of compatible pens is provided in these instructions for use.
- Patients with mental illness, convulsions/seizures or other physical impairments and children might need to be assisted while using the device.
- Keep out of reach of small children as the product contains small parts that could be swallowed.
- KD-*PENOFINE* is non-toxic, sterile and non-pyrogenic within shelf life if the individual package has not been opened or damaged.
- Do not use after the expiry date. If used after the expiry date the device may not be sterile.
- Visually inspect and carefully check the product and individual package integrity and also the expiry date before use. Improper transport and handling may cause structural and/or functional damage to device or packaging.
- Do not use if the product or individual package is impaired, damaged, opened or contaminated as this can alter the device's sterile state.
- Use the product immediately after opening the individual package.
- Do not bend the needle. If the needle is bent or damaged, do not try to straighten or use the product as it may lead to needlestick injuries / impaired integrity of the needle.
- Handle and dispose of with care to avoid needlestick injuries.
- The pen needle is for single use on a single patient only — always use a new needle for each injection and discard it directly afterwards. Re-use of this single-use device might lead to infections as a used needle is no longer sterile. Furthermore, injections with an already used needle might be more painful as the lubricant is removed and the tip shape is deformed.
- Do not clean or re-sterilize the product.
- Never store the insulin pen with a mounted pen needle as the mounted needle provides an open channel to the insulin medication inside the cartridge - the insulin medication may be altered (e.g. crystallization) or contaminated, air may entry into the cartridge and dosage might become inaccurate.

Instructions for use:

Before using the pen needle, consult and strictly follow the insulin pen's instructions for use. Always use the needle size, injection technique and injection site recommended by your doctor for administration.

1. Open the sealing strip of the individual package (see Fig.1);
2. Attach the product to a compatible insulin pen by screwing it clockwise onto the pen (see Fig.2);
3. Remove both, the outer big protective cap (see Fig.3) and the inner small protective cap (see Fig.4), before an injection. Do not dispose of the big protective cap so you can replace it in step 7.
Warning: If both, the big protective cap and the small protective cap, are not removed before use, the medication or dose might not be injected, which can result in harm.
4. Perform a priming test prior to use to remove any air and to ensure correct dose administration according to the insulin pen's instructions for use (see Fig.5). Always keep backup pen needles as a replacement. Make sure that the injection dosage is set correctly (see Fig.6).
5. Disinfect the skin of the injection site.
6. Hold the insulin pen in a 90° angle to the skin surface at the insertion site and insert the needle vertically. Once the pen needle is correctly inserted, press the injection button all the way down, being careful not to change the angle of the insulin pen (see Fig.7). Keep pressing down the injection button. After administration keep the needle inserted for at least 10 seconds before removing the needle.
7. When the injection has been completed, replace the big protective cap to cover the needle (see Fig.8) and rotate the pen needle anticlockwise to detach it from the insulin pen (see Fig.9).
Note: If after the needle is removed from the skin you notice medication dripping, keep the needle inserted longer the next time, at the same dosage.
8. Dispose of the used pen needle in accordance with applicable regulations and guidelines.

Disposal:

Dispose of the used pen needle after single use following local regulations and guidelines for sharps disposal by using an appropriate container for sharp medical waste.

Storage:

Store in a cool and dry place. Do not expose to heat or direct sunlight.

Notice:

Any serious incident that has occurred in relation to the device must be reported to KD Medical GmbH Hospital Products and to the competent authority of the EU Member State (or state with the same rules (Regulation (EU) 2017/745 on medical devices)) in which the user and/or patient is established.

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