# KD-PENOFINE S

## Safety Pen Needles

Sicherheits-Pen-Nadeln • Säkerhetspennkanyler • Sikkerhedspennåle • Sikkerhetspennkanyler • Turvakynäneulat • Veiligheidspennaalden • Aiguilles pour stylos avec dispositif de sécurité • Bezpieczne igły do pena • Bezpečnostní jehly pro injekční pera • Bezpečnostné ihly na perá • Sigurnosne igle za inzulinsku olovku • Ace stilou injector cu mecanism de siguranță • Безопасни игли за писалки

## Compatible pens






Kompatible Pens • Kompatibla pennor • Kompatible penne • Kompatible penner • Yhteensopivat kynät • Compatibele pennen • Stylos compatibles • Kompatybilne peny • Kompatibilní pera • Kompatibilné perá • Kompatibilne inzulinske olovke • Ace stilou injector compatibile • Съвместими писалки

|                |                                                                                                                                                          |
|----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|
| BD Medical     | Puregon Pen®                                                                                                                                             |
| Berlin Chemie  | BerliPen® areo 2,<br>BerliPen® areo 3                                                                                                                    |
| Bioton         | GensuPen Improve,<br>GensuPen2                                                                                                                           |
| Eli Lilly      | Humalog® KwikPen™,<br>HumaPen®, HumaPen® Ergo,<br>HumaPen® Ergo II, HumaPen® Luxura,<br>HumaPen® Luxura HD, Saline KwikPen™,<br>HumaPen® Savvio          |
| Gan & Lee      | GanleePen™                                                                                                                                               |
| Novo Nordisk   | NovoMix®30 FlexPen®,<br>NovoPen® 3, NovoPen® 4, NovoPen® 5,<br>NovoPen® 6, NovoPen® Junior,<br>NovoPen Echo®, NovoPen Echo® Plus,<br>Tresiba® FlexTouch™ |
| Owen Mumford   | AutoPen®, AutoPen® 2,<br>AutoPen® 24                                                                                                                     |
| Sanofi Aventis | ALLSTAR®, ALLSTAR® Pro,<br>JuniorSTAR®, Lantus® SoloStar®,<br>OptiPen® Pro 1, Sanofi Aventis Lantus®,<br>TactiPen®                                       |
| Ypsomed        | YpsoPen®                                                                                                                                                 |

EN ISO 11608-2

## Figures 1-6

Abbildung 1 bis 6 • Figur 1-6 • Figur 1-6 • Figur 1-6 • Kuvat 1-6 • Afbeelding 1-6 • Figures 1-6 • Rysunki 1-6 • Obrázky 1-6 • Obrázky 1 – 6 • Slike 1-6 • Figurile 1-6 • Фигури 1 – 6

# en Safety Pen Needles

## Brief description:

KD-*PENOFINE S* pen needles consist of the following components: outer sheath, sealing strip, shielding protective cap, sliding sleeve, spring and the actual needle fixed in a plastic hub.

The safety mechanism of KD-*PENOFINE S* covers the needle after use. In the activated position, the safety mechanism guards against accidental needlestick injuries during handling of the used product and its disposal.

## Materials used:

Stainless steel, polypropylene (PP), ABS, silicone oil (lubricant), UV glue (adhesive)

## Intended purpose:

KD-*PENOFINE S* pen needles are single use medical devices intended to be connected to a compatible insulin pen for subcutaneous insulin administration within the frame of a diabetes mellitus treatment regimen. KD-*PENOFINE S* pen needles are sterilized using ethylene oxide.

## Indication:

Subcutaneous injections of insulin with a compatible insulin pen

## Contraindications:

- Insulin/insulin derivatives injections during episodes of hypoglycemia
- Hypersensitivity to insulin/insulin derivatives or one of its excipients
- Inflammation, tumors, ulceration or other pathologic skin / subcutaneous tissue alterations or lesions such as bruising or indurations or scar tissue at the injection site

## Residual risks, undesirable side effects:

- Injection-related side effects such as bleeding, bruising (hematoma) or skin reactions (e.g. redness)
- Lipohypertrophy
- Hypoglycemia
- Needle breaking or bending
- Pain
- Insulin backflow (leakage)

## Intended user:

The use of this product is restricted to qualified healthcare professionals and trained laypersons such as patients or caregivers.

## Intended patient population:

KD-*PENOFINE S* pen needles are intended for patients who need insulin injection for blood glucose control without gender or age constraints. The different needle sizes available are to be selected by experienced healthcare professionals regarding suitability for the intended injection sites and the patient's body condition.

## ⚠ Warnings, precautions:

- Read instructions before use. KD-*PENOFINE S* is to be used according to these instructions for use.
- **As a layperson, consult your doctor regarding the most suitable needle size, injection technique and injection sites to be used and always follow your doctor's instructions.**
- Check the compatibility between the pen needle and the insulin pen — a list of compatible pens is provided in these instructions for use.
- Patients with mental illness, convulsions/seizures or other physical impairments and children might need to be assisted while using the device.
- Keep away from small children as the product contains small parts that could be swallowed.
- KD-*PENOFINE S* is non-toxic, sterile and non-pyrogenic within shelf life if the individual package has not been opened or damaged.
- Do not use after the expiry date. If used after the expiry date the device may not be sterile.
- Visually inspect and carefully check the product and individual package integrity and also the expiry date before use. Improper transport and handling may cause structural and/or functional damage to device or packaging.
- Do not use if the product or individual package is impaired, damaged, opened or contaminated as the device is no longer sterile in this case.
- Use the product immediately after opening the individual package.
- Do not bend the needle. If the needle is bent or damaged, do not try to straighten or use the product as it might lead to needlestick injuries / impaired integrity of the needle.
- Handle and dispose of with care to avoid needlestick injuries.
- The pen needle is for single use on a single patient only — always use a new needle for each injection and discard it directly afterwards. Re-use of this single-use device might lead to infections as a used needle is no longer sterile. Furthermore, injections with an already used needle might be more painful as the lubricant is removed and the tip shape is deformed.
- Do not clean or re-sterilize the product.
- Never store the insulin pen with a mounted pen needle as the mounted needle provides an open channel to the insulin medication inside the cartridge - the insulin medication may be altered (e.g. crystallization) or contaminated, air may entry into the cartridge and dosage might become inaccurate.

## Instructions for use:

Before using the pen needle, consult and strictly follow the insulin pen's instructions for use. Always use the needle size, injection technique and injection site recommended by your doctor for administration.

1. Open the sealing strip of the individual package (see Fig.1);
2. Attach the product to a compatible insulin pen by screwing it clockwise onto the pen (see Fig.2);
3. Remove the outer sheath (see Fig. 3).
4. Perform a priming test prior to use to remove any air and to ensure correct dose administration. Make sure that the injection dosage is set correctly.
5. Disinfect the skin of the injection site.
6. Hold the insulin pen in a 90° angle to the skin surface at the insertion site and insert the needle vertically. Once the pen needle is correctly inserted, press the injection button all the way down, being careful not to change the angle of the insulin pen. Keep pressing down the injection button. After administration keep the needle inserted for at least 10 seconds before removing the needle (see Fig.4).
7. When the injection has been completed, remove the insulin pen. The safety protective mechanism will activate automatically to shield the needle point (see Fig.5): upon hearing a clear "click" sound, the needle point is shielded, and the color indicator in the visual window of the shielding protective cap disappears:

Color indicator is visible in the visual window:

Safety mechanism is not activated.

Color indicator has disappeared: Safety mechanism is activated.

Note: If after the needle is removed from the skin you notice medication dripping, keep the needle inserted longer the next time, at the same dosage.

8. Rotate the pen needle anticlockwise to detach it from the Insulin pen (see Fig.6). Dispose of the used pen needle in accordance with applicable regulations and guidelines.

## Disposal:

Dispose of after single use following applicable regulations and guidelines for sharps disposal by using an appropriate container for sharp medical waste.

## Storage:

Store in a cool and dry place. Do not expose to heat or direct sunlight.

## Notice:

Any serious incident that has occurred in relation to the device must be reported to KD Medical GmbH Hospital Products and to the competent authority of the EU Member State (or state with the same rules (Regulation (EU) 2017/745 on medical devices)) in which the user and/or patient is established.

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